


#### Resona 7

Premium Ultrasound System

New Waves in Ultrasound Innovation







# It rises.

### With core platform advantages of ZST<sup>+</sup>

The channel data based ZST<sup>+</sup> is an extraordinary innovation, representing an ultrasound evolution. Transforming ultrasound metrics from conventional beamforming to channel data based processing, ZST<sup>+</sup> is able to deliver multiple imaging advances: Advanced Acoustic Acquisition, Dynamic Pixel Focusing, Sound Speed Compensation, Enhanced Channel Data Processing and Total Recall Imaging.





#### **Advanced Acoustic Acquisition**

By transmitting and receiving a relatively smaller number of large zones, Advanced Acoustic Acquisition extracts more information from each acquisition, 10 times faster than a conventional line-by-line beamforming method.

# Line-by-line Acquisition Advanced Acoustic Acquisition Time to Form Ultrasound Image Additional Time for Advanced Modes


#### **Dynamic Pixel Focusing**

Dynamic Pixel Focusing technology allows the Resona 7 to achieve extreme uniformity in pixel level throughout the whole field of view. Now there's no need to adjust the focal positions to achieve uniformity across patient exams.




Actual Tx Profile

Dynamic Pixel Focusing Based On ZONE Sonography®



**Sound Speed Compensation** 

By retrospectively analyzing complete channel data stored in channel data memory, the Resona 7 is able to intelligently choose the optimal sound speed to improve image accuracy even with tissue variation, allowing for adaptive tissue-specific optimization.



SSC OFF

Effective Tx Profile

**Enhanced Channel Data Processing** 

Channel data based ZST<sup>+</sup> provides Enhanced Channel Data Processing for greatly improved imaging clarity. By multiple and retrospective channel data processing, it makes the best use of acoustic information for image improvement.

- HD Scope: higher definition image within ROI.
- Coherent Spatial Synthesis: further improved image quality of spatial compound.








SSC ON

HD Scope ON

#### **Total Recall Imaging**

As ZST<sup>+</sup> captures and stores the complete acoustic raw data set, Total Recall Imaging allows system to do retrospective processing on channel data and also permits users to modify numerous imaging parameters on stored images to maximize clinical output.





Testicle Perfusion




Ankle Trauma



Follicles




CEUS of Liver Lymphoma



Elastography of Thyroid Mass




CCA and Jugular Vein



3D Uterus Septus




3D Fetal Face



**Aortic Regurgitation** 



3D Fetal Spine



TT QA

# It releases.

## A new standard of image clarity


Better vision, deeper understanding. Based on the cutting-edge ZST<sup>+</sup> platform, Resona 7 redefines a new standard of image performance to meet the needs of the most challenging clinical practices.

# It progresses.


# Innovative clinical tools for confident diagnosis

#### **V Flow**

V Flow (Vector Flow) is a novel approach for vascular hemodynamic analysis. V Flow uses color coded vector arrows to indicate the velocity's magnitude and direction of blood cells. With over 300 frames per second, it provides extremely vivid, accurate and angle-independent visualization of complex vascular hemodynamics profiles. With comprehensive data information, V Flow is the most valuable tool for vascular clinical research.




V Flow of Carotid Bulb and JV




V Flow of CCA and ICA

#### **UWN<sup>+</sup> Contrast Imaging**

UWN<sup>+</sup> (Ultra-Wideband Non-linear Plus) CEUS enables the Resona 7 to detect and utilize both 2<sup>nd</sup> harmonic and non-linear fundamental signals, generating significantly enhanced images, resulting in greater sensitivity of minor signals and longer agent duration with lower MI.



TIC Analysis of CEUS




Micro Flow Enhancement




#### **iFusion with Respiration Compensation**

Bringing the precision of fusion imaging to a new level, Mindray's pioneering, innovative and exclusive respiration compensation technology - supported by a sensitive magnetic motion sensor with millimeter accuracy - can help eliminate distortion and fusion inaccuracy caused by patient respiration.



iFusion without Respiration Compensation



iFusion with Respiration Compensation



# It leads.

# Forwarding smart to clinical intelligence

The Resona 7 elevates clinical intelligence to a new level with a complete solution that enables clinicians to manage both routine and advanced studies more efficiently, consistently, and accurately, from acquisition to calculation. As an example, Smart Planes shows exceptional intelligence in accurate diagnosis and analysis of fetal central nervous system (CNS).

#### **Smart Planes**

Mindray's exclusive pioneering technology positions the Resona 7 as the industry's first ultrasound system to allow fully automatic and accurate detection of the most significant planes and frequently used measurements of fetal CNS, leading to intelligent diagnosis, improved throughput, and reduced user dependency.

Smart Planes provides a user-friendly tool that greatly improves scanning efficiency through increased accuracy coupled with automated operation. With a simple button click on a 3D fetal brain volume image, the standard CNS scanning planes (MSP, TCP, TTP and TVP) and a range of related anatomical measurements (BPD, HC, OFD, TCD, CM and LVW) are obtained immediately.

Dandy Walker & ACC






Smart Planes

MRI

Abnormal CM with Suspected Dandy Walker Syndrome





Smart Planes


MRI

#### **Smart FLC**

Smart FLC automatically detects the number of follicles and calculates each volume from a 3D ovarian volume image, assuring accurate assessment of follicles, especially with IVF exams.

#### **Smart OB/NT**

Automatic measurements of the most frequently examined parameters, including BPD, HC, FL, AC, OFD and even NT as early as 11 weeks, are available with a single click for higher productivity and reproducibility.



Smart FLC



Smart N

# It senses.

### Ensuring a better user experience

The Resona 7 is designed around you. Gesture-based operation opens up a new trend in cart-based ultrasound with an agile, smart, and intuitive user experience beyond your expectations. A six-direction floating control panel with electronic height adjustment provides scanning comfort in any position. Inspired innovations drive a better user experience.











\_\_\_\_ 21.5"
high resolution LED monitor

tilting multi-gesture touch screen


—— 6
direction floating control panel

\_\_\_\_ Gel warmer with temperature control

\_\_\_\_ Pinless transducer with light indicator

— Central and swivel lock







With over 20 years of experience, Mindray hosts a wide range of ultrasound imaging solutions including cart-based and portable systems. Being exported to over 190 countries, Mindray ultrasound systems are today being used by medical professionals for general as well as highly dedicated clinical utility. With a global R&D base spanning over Asia, Europe and America, the ultrasound solutions by Mindray are a result of an integral cooperation with the medical community, allowing for the ultrasound systems to be extremely user centric in terms of performance and usability. Mindray is well positioned to become one of the leading ultrasound imaging solutions provider.

Mindray Building, Keji 12th Road South, High-tech Industrial Park, Nanshan, Shenzhen 518057, P.R. China Tel: +86 755 8188 8998 Fax: +86 755 26582680 E-mail: intl-market@mindray.com www.mindray.com

Mindray is listed on the NYSE under the symbol "MR"  $\,$ 

mindray | Numbers with made are registered trademarks or trademarks owned by Shenzhen Mindray Bio-medical Electronics Co., LTD.

© 2015 Shenzhen Mindray Bio-Medical Electronics Co., Ltd. All rights reserved. Specifications subject to changes without prior notice.

P/N:ENG-Resona7-210285x16P-20151105

